

CHEMISTRY 203

CHEMICAL EQUILIBRIUM

Determination of an Equilibrium Constant


I. Purpose:

1. To review the concepts and principles of **chemical equilibrium**.
2. To study the **effects of temperature and concentration changes** on the position of equilibrium of a reversible reaction.
3. To learn the principles and use of a **spectrophotometer**.
4. To apply the principles of chemical analysis by **spectrophotometric** means.
5. To determine, spectrophotometrically, the **equilibrium constant K_c** for the reaction between iron (III) ion and thiocyanate ion (SCN^-) at a given temperature.

II. Theory:

A. Spectrophotometry

- A *spectrophotometric method* of analysis involves the interaction of *electromagnetic radiation* (EM) with matter.
- The most common regions of EM spectrum used for analysis are the **ultraviolet UV**, **visible**, and the **infrared IR** regions.
- When a substance is colored it absorbs light in the visible region.
- **Visible region** has a wavelength ranging from **400-700 nm**.



a. What do we see?

- When continuous **electromagnetic radiation** (light) passes through a substance, a portion of the light may be absorbed. The remaining light exits the sample and when passed through a prism will yield a spectrum with gaps in it, this is called the *absorption spectrum*.
- The color the **eye sees** is due to the wavelengths of light the sample did not absorb, that is we see the **transmitted** color.

b. Colors and wavelengths in the visible region of the electromagnetic spectrum:

Color absorbed	Wavelength (nm)	Color transmitted
Red	750-610	Green-blue
Orange	610-595	Blue-green
Yellow	595-580	Violet
Green	580-500	Red-violet
Blue	500-435	Orange yellow
Violet	435-380	yellow

c. Relationship between energy, frequency and wavelength:

$$E = h\nu \qquad \nu = c / \lambda$$

h = Planck's constant = 6.626×10^{-34} Js

C = speed of light = 3.00×10^8 m s⁻¹

λ = wavelength in m, 1 nm = 10^{-9} m

ν = frequency s⁻¹ or hz

E decreases ↓ as λ increases ↑

d. Absorption Spectrum of a Chemical Substance and Spectrophotometer

Absorption spectrum: A plot of absorbance **A** versus wavelength **λ**

- The wavelength (**λ_{max}**) at which the **maximum** absorption of the (EM) radiation by the absorbing substance occurs is **determined** and **set** on the **spectrophotometer**.
- **Spectrophotometer:** A laboratory instrument that measures the amount of light transmitted through a sample.



e. Spectrophotometer Components




Figure 1

- **Source of light**, stable source of radiant energy
- **Wavelength selector** or **monochromator**, that isolates a limited region of the spectrum for measurement
- **Sample container**.
- **Radiation detector**, which converts radiant energy to a measurable electrical signal.
- **Recorder**, signal processing and read out unit.

f. Transmittance and Absorbance

The **transmitted** light $I_t < I_0$ **incident** light

Transmittance: $T = \frac{I_t}{I_0} < 1$

$$\% T = \frac{I_t}{I_0} \times 100$$

Absorbance: $A = -\log T$

$$A = -\log \frac{\% T}{100}$$

g. Beer-Lambert's Law:

- The wavelength with which the compound interacts with light is an experimental parameter.
- Several factors control the amount EM radiation (light energy) that a sample absorbs.
 1. **Concentration** of the absorbing substance.
 2. **Thickness** of the sample containing the absorbing substance (determined by the width of the cuvet).
 3. Probability of light absorption by the absorbing substance (called the **molar absorptivity coefficient** or the **extinction coefficient**).


- A spectrophotometer measures the amount of light absorbed by a colored substance. **Beer-Lambert's Law** relates quantitatively the extent of absorption to the concentration of the absorbing species according to the following relationship:

$$A = \epsilon \cdot b \cdot c$$

- A: optical activity or absorbance at specific wavelength λ_{\max}
- ϵ : extinction coefficient or molar absorptivity ($M^{-1}cm^{-1}$), f (λ , temp. solvent)
- b: width of cuvet (cell), light path length, or depth(cm)
- c: Molar concentration (M)

Beer-Lambert's law applies best near maximum absorption regions, and for dilute solutions.

A plot of absorbance versus concentration data is linear; it provides a **calibration curve**. The calculated slope and the absorbance data of an absorbing species can be used to determine the molar concentration of this species in an unknown solution. $\text{slope} = \epsilon \cdot b \text{ (M}^{-1}\text{)}$



pharmlabs.unc.edu/spectrophotometry/text.htm

II. Theory (continued):

B. Chemical Equilibrium

- Few chemical reactions proceed in only one direction. Most are **reversible**.
- A reversible chemical reaction leads to an **equilibrium state** which is characterized by **constant composition**.
- Although the reaction has not stopped at equilibrium, the **rate of disappearance of the reactants is equal to the rate of formation of reactants from the products**.
- Chemical equilibrium is a **dynamic** process.
- Chemical equilibrium represents a **balance** between **forward** and **reverse** reactions.
- Thus, the **concentrations of reactants and products remain constant**.

• **Homogeneous** equilibrium applies to reactions in which all reacting species are in the **same phase**.

• In general, for a reversible reaction at a specific temperature,

$$aA + bB \rightleftharpoons cC + dD$$

• The molar concentrations (mol/L) of reactants and products at equilibrium are related by the **equilibrium constant expression K_c** :

$$K_c = \frac{[C]^c [D]^d}{[A]^a [B]^b}$$

Note: **The equilibrium mixtures can vary, but *the value for K at a specific temperature is constant.***

Le Chatelier's Principle


Le Chatelier's principle summarizes the effect of stressing a chemical system at equilibrium as follows:

When a **stress** is applied to a system at equilibrium the system will **shift** in such a direction to **relieve** the stress and establish a **new** state of equilibrium.

Factors that affect chemical equilibrium:

- Change in concentration
- Change in temperature
- Change in volume and pressure

Formation of FeSCN²⁺ complex ion:



The equilibrium constant is defined as:

$$K_c = \frac{[\text{FeSCN}^{2+}]}{[\text{Fe}^{3+}][\text{SCN}^{-}]}$$

FeSCN²⁺ complex ion has a deep blood red color and absorbs light in the visible region at $\lambda_{\text{max}} = 447 \text{ nm}$


In this experiment the concentration of FeSCN²⁺ will be determined spectrophotometrically, using Beer-Lambert's law.

$$A = \epsilon b c$$

- **Part One.** A set of standard solutions will be prepared, their absorbance measured, and the absorbance for each solution is plotted vs. the molar concentrations of FeSCN²⁺ to establish the calibration curve.
- In preparing the standard solution, the Fe³⁺ concentration will far exceed the SCN⁻ concentration.
- This huge excess of Fe³⁺ pushes the equilibrium far to the right consuming all the SCN⁻ placed in the system .
- As a result, the FeSCN²⁺ concentration at equilibrium approximates the original SCN⁻ concentration.
- We assume that the position of the equilibrium is driven so far to the right that all the SCN⁻ is complexed forming FeSCN²⁺.

- **Part Two.** The chemical system is prepared by mixing known molar concentrations of Fe^{3+} and SCN^- , so by knowing the initial concentrations and measuring the equilibrium concentration of FeSCN^{2+} spectrophotometrically, the equilibrium concentrations of Fe^{3+} and SCN^- are calculated.

III. Calculations:



At equil. $(n_{\text{o Fe}^{3+}} - n_{\text{FeSCN}^{2+}})$ $(n_{\text{o SCN}^-} - n_{\text{FeSCN}^{2+}})$ $n_{\text{FeSCN}^{2+}}$

Where: $n_{\text{o Fe}^{3+}} = (M_{\text{o Fe}^{3+}}) (V_{\text{o Fe}^{3+}})$

$n_{\text{o SCN}^-} = (M_{\text{o SCN}^-}) (V_{\text{o SCN}^-})$

$n_{\text{FeSCN}^{2+}} = x = [\text{FeSCN}^{2+}]_{\text{eq}} V_{\text{tot}} = (A/\epsilon b) V_{\text{tot}}$

A is the measured absorbance for that particular solution.

ϵb = slope of calibration curve

$$n_{\text{eq Fe}^{3+}} = n_{\text{o Fe}^{3+}} - n_{\text{FeSCN}^{2+}}$$

$$n_{\text{eq SCN}^{-}} = n_{\text{o SCN}^{-}} - n_{\text{FeSCN}^{2+}}$$

$$[\text{Fe}^{3+}]_{\text{eq}} = n_{\text{eq Fe}^{3+}} / V_{\text{tot}}$$

$$[\text{SCN}^{-}]_{\text{eq}} = n_{\text{eq SCN}^{-}} / V_{\text{tot}}$$

$$K_c = \frac{[\text{FeSCN}^{2+}]}{[\text{Fe}^{3+}][\text{SCN}^{-}]}$$