

29.5 / 30



Chemistry 205 Report

Redox Titration

Name:

Date: 22-02-2013

**Purpose:** 1. To learn a technique in volumetric analysis Redox titration  
 2. To review the stoichiometry of an oxidation-reduction reaction  
 3. To determine the concentration of an unknown sodium oxalate ( $\text{Na}_2\text{C}_2\text{O}_4$ ) solution by titrating it against standardized potassium permanganate solution ( $\text{KMnO}_4$ )  
 4. To determine the percent by mass of Fe (II) in the form of ferrous ammonium sulfate  $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$  in a mixture by redox titration.

**Part I**  
**Table 1. Preparation of Standard  $\text{Na}_2\text{C}_2\text{O}_4$  Solution**

Mass of bottle + cap + $\text{Na}_2\text{C}_2\text{O}_4$	10,0868 g $\pm 0.0001$ g
Mass of empty bottle + cap	9,2768 g
Mass of $\text{Na}_2\text{C}_2\text{O}_4$	0,8100 g
Volume of solution	250,0 ml $\pm 0.1$ ml

$$M_{\text{Na}_2\text{C}_2\text{O}_4} = \frac{m}{V} = \frac{0,8100 \text{ g}}{134,00 \text{ ml}} = 6,045 \times 10^{-3} \text{ M}$$

$$M_{\text{Na}_2\text{C}_2\text{O}_4} = \frac{m_{\text{Na}_2\text{C}_2\text{O}_4}}{V_{\text{Na}_2\text{C}_2\text{O}_4}} = \frac{0,8100 \text{ g}}{250,0 \text{ ml}} = 0,00324 \text{ M}$$

**Table 2. Standardization of  $\text{KMnO}_4$  Solution**

Trial	Buret upper reading	Buret lower reading	Volume of $\text{KMnO}_4$ $\pm 0,02$ ml
1	0,00	26,60 $\pm 0,2$	26,60
2	0,00	26,80	26,80
3	0,00	26,20	26,20
Blank	0,00	0,40	0,40

$$\bar{V}_{\text{KMnO}_4} = 26,53$$

$$V_{\text{KMnO}_4} = \bar{V}_{\text{KMnO}_4} - V_{\text{Blank}} = 26,53 - 0,4 = 26,13 \text{ ml}$$

Ionic equation:



Volume of  $\text{Na}_2\text{C}_2\text{O}_4 = 25,00 \text{ ml}$ .

Calculations: (show your work)

$$M \text{ of } \text{KMnO}_4 = 9,254 \times 10^{-3} \text{ M}$$

$$\frac{n_{\text{MnO}_4^-}}{2} = \frac{n_{\text{C}_2\text{O}_4^{2-}}}{5}$$

$$\frac{(M \times V)_{\text{MnO}_4^-}}{2} = \frac{(M \times V)_{\text{C}_2\text{O}_4^{2-}}}{5}$$

$$M_{\text{Na}_2\text{C}_2\text{O}_4} = 2 \times 0,02418 \times 25,0 \times 10^{-3} = 9,254 \times 10^{-3} \text{ M}$$





**Table 3.** Analysis of Unknown Oxalate

Unknown #: 5

Trial	Buret upper reading	Buret lower reading	Volume of KMnO <sub>4</sub> ml
1	0,00	27,50	27,50
2	0,00	27,30	27,30
3	0,00	27,10	27,10

$$\bar{V}_{\text{KMnO}_4} = 27,30$$

$$V_{\text{KMnO}_4} = \bar{V}_{\text{KMnO}_4} - V_{\text{Blank}}$$

$$= 27,30 - 0,40$$

$$= 26,9 \text{ ml.}$$

Volume of unknown Na<sub>2</sub>C<sub>2</sub>O<sub>4</sub> = 25.00 ml ± 0.02 ml

**Calculations:** (show your work)

Molarity of unknown Na<sub>2</sub>C<sub>2</sub>O<sub>4</sub> =

$$\frac{(M \times V)_{\text{H}_2\text{O}_4^{2-}}}{2} = \frac{(M \times V)_{\text{Na}_2\text{C}_2\text{O}_4}}{5}$$

$$M_{\text{Na}_2\text{C}_2\text{O}_4} = \frac{5 \times 9,254 \times 10^{-3} \times 26,9 \times 10^{-3}}{2 \times 25,00 \times 10^{-3}} = 0,02489 \text{ M}$$



% error = 2% ✓

**Part II**

**Table 4.** Preparation of the Iron Unknown Solution

Unknown #: 12

Mass of bottle + cap + impure Fe(NH <sub>4</sub> ) <sub>2</sub> (SO <sub>4</sub> ) <sub>2</sub> ·6H <sub>2</sub> O	12,3927 ± 0,001
Mass of empty bottle + cap	10,3689
Mass of impure Fe(NH <sub>4</sub> ) <sub>2</sub> (SO <sub>4</sub> ) <sub>2</sub> ·6H <sub>2</sub> O	2,0238g
Volume of solution	100.0 ml ± 0.1 ml

**Table 5.** Titration of the Iron Unknown Sample with KMnO<sub>4</sub> Solution

Trial	Buret upper reading	Buret lower reading	Volume of KMnO <sub>4</sub>
1	0,00	5,80	5,80
2	5,80	11,70	5,90
3	11,70	17,30	5,60

$$\bar{V}_{\text{KMnO}_4} = 5,77 \text{ ml}$$



Ionic equation:



Calculations: (show your work)

Number of moles of  $\text{Fe}^{2+}$  (in titrated sample) =  $\frac{n_{\text{Fe}^{2+}}}{5} = \frac{n_{\text{MnO}_4^-}}{1}$

$$5(\text{H} \times \text{V})_{\text{MnO}_4^-} = \frac{1}{5} \times 9,254 \times 10^{-3} \times 5,77 \times 10^{-3}$$

$$= 2,67 \times 10^{-4} \text{ mol.}$$

Number of moles of  $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$  (in titrated sample) =  $n_{\text{Fe}^{2+}} = 2,67 \times 10^{-4} \text{ mol.}$

Number of moles of  $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$  (in the original sample) =  $n_{\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}} \times 10$

Mass of  $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O} = n \times M = 2,67 \times 10^{-3} \times 395,00 = 2,67 \times 10^{-3} \text{ mol.}$

$$= 1,05 \text{ g}$$

% by mass of  $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O} = \frac{\text{mass of iron}}{\text{mass of powder}} \times 100 = \frac{1,05}{2,0238} \times 100$

good job

% error = 2,8%

$$= 51,9\%$$

Questions

1. Give the oxidation number of the underlined atoms in the following species:

- a.  $\underline{\text{C}}_2\text{O}_4^{2-}$   $2\text{C} + 4(-2) = -2 \Rightarrow \text{C} = +3$
- b.  $\text{Mg}_3\underline{\text{N}}_2$   $3(+2) + 2\text{N} = 0 \Rightarrow \text{N} = -3$
- c.  $\text{H}_3\underline{\text{As}}\text{O}_3$   $3(+1) + \text{As} + 3(-2) = 0 \Rightarrow \text{As} = +3$
- d.  $\underline{\text{W}}\text{O}_4^{2-}$   $\text{W} + 4(-2) = -2 \Rightarrow \text{W} = +6$
- e.  $\underline{\text{Pt}}\text{Cl}_6^{2-}$   $\text{Pt} + 6(-1) = -2 \Rightarrow \text{Pt} = +4$
- f.  $\underline{\text{Sb}}\text{F}_6^-$   $\text{Sb} + 6(-1) = -1 \Rightarrow \text{Sb} = +5$



2.  $\text{SO}_2$  present in air is mainly responsible for the acid rain phenomenon. Its concentration can be determined by titrating against a standard permanganate solution in acidic media as follows:



- a. Balance the above reaction, show your work.
- b. Calculate the number of grams of  $\text{SO}_2$  in a sample of air if 7.37 ml of 0.00800 M  $\text{KMnO}_4$  solution are required for titration.

