

Chemistry 101 Laboratory
Fall 2005 - 2006

Lecture 5
Acid - Base Titration

purpose


- To learn the **concept** and **technique** of titration.
- To **standardize** a sodium hydroxide (NaOH) solution against a primary standard acid.
- To determine the **concentration** of an **unknown acid** by titration with the standardized base solution.

Titration

In **titration** a solution of accurately known concentration, **called a standard solution**, is added gradually to another solution of unknown concentration (or vice versa) until the chemical reaction between the two solutions is **complete**.

Equivalence point – the point at which the reaction is **complete**.

Indicator – substance that changes color at (or near) the equivalence point



Slowly add base
to unknown acid
UNTIL

the indicator
changes color







Titration (Cont'd)

- ***End Point:*** the point in the titration at which the indicator changes color (visually determined).
- ***Conditions for a good titration***
 - Rapid and complete reaction.
 - Reaction of known stoichiometry (no side products).
 - End point easily detected (use proper indicator).

Acid – Base Titration

Involves a **neutralization reaction** which is the complete reaction between an acid and a base.





What volume of a 1.420 M NaOH solution is Required to titrate 25.00 mL of a 4.50 M H₂SO₄ solution?



WRITE THE CHEMICAL EQUATION!



$$n_{\text{NaOH}} = 2 n_{\text{H}_2\text{SO}_4}$$

$$(M \times V)_{\text{NaOH}} = 2 (M \times V)_{\text{H}_2\text{SO}_4}$$


$$1.420 \text{ mol/L} \times V_{\text{NaOH}} = 2 (4.50 \text{ mol/L} \times 25.00 \times 10^{-3}\text{L})$$

$$V_{\text{NaOH}} = 0.158 \text{ L} = 158 \text{ mL}$$

Standard Solution

- Is a solution of ***accurately known concentration***.
- prepared by dissolving an exact amount of the solute, followed by dilution, to form a definite volume of solution.

Preparing a Solution of Known Molarity




Standard solution (cont'd)


- ***Properties of a good primary standard***
 - High molar mass
 - Stable
 - Not hygroscopic
 - Highly pure, cheap and available
- ***Example:*** potassium hydrogen phthalate,
 $\text{KHC}_8\text{H}_4\text{O}_4$, abbreviated as KHP,
molar mass = 204.23 g / mol.

Standard solution (cont'd)


- Sodium hydroxide is not a good primary standard solution. Solid NaOH absorbs water from air (hygroscopic), and its solution reacts with carbon dioxide.
- In this experiment a sodium hydroxide (NaOH) solution will be **standardized** by titration with a **primary standard acid (KHP)** as follows:



or



or



net ionic equation

Experimental

Three steps are involved:

- 1- Preparation of a primary standard potassium hydrogen phthalate (KHP) solution.
- 2- Standardization of 0.1 M sodium hydroxide solution with KHP.
- 3- Titration of an unknown monoprotic acid with the standardized sodium hydroxide solution.

1- Preparation of KHP Primary Standard

- Weigh KHP bottle using the analytical balance.
- Transfer the KHP to a 250 mL volumetric flask using a funnel.
- Add some water, dissolve, dilute and shake, make up to the mark and homogenize.
- Weigh the empty bottle
- Calculate molarity of KHP

$$M = n / V = m / (\text{molar mass}) \times V$$

$$M = m \text{ (g)} / 204.22 \text{ g/mol} \times (250.0 \times 10^{-3} \text{ L})$$

2- Standardization of 0.1M NaOH

- Rinse the buret with NaOH.
- Fill it to the mark, make sure that there are no air bubbles in the tip.
- Pipet 10 mL of KHP into an Erlenmeyer flask and add two drops of phenolphthalein indicator.
- Add 10mL of water.
- Titrate against KHP to the end point (light pink color).
- Repeat the titration 3 times.
- Calculate M of NaOH.

$$(M \times V) \text{ of NaOH} = (M \times V) \text{ of KHP}$$

3- Titration of an Unknown Acid

- Proceed as in part 2.
- Use the unknown instead of KHP
- Titrate against the standardized NaOH solution using phenolphthalein as an indicator.
- Repeat three times.
- Calculate the molarity of the unknown.