

# Chemistry 101 Laboratory

## Fall 2005 - 2006


### Lecture 9

### Solubility as a Function of Temperature

# Purpose

- To determine the **solubility** of oxalic acid ( $\text{H}_2\text{C}_2\text{O}_4$ ) in water at three different temperatures by titrating **saturated solutions** of the acid with a ***standard*** solution of NaOH.

## Equation Involved and Calculations



Moles of  $\text{H}_2\text{C}_2\text{O}_4 = \text{moles of NaOH} / 2$

Moles of  $\text{H}_2\text{C}_2\text{O}_4 = (\text{M} \times \text{V}) \text{ NaOH} / 2$

Mass of  $\text{H}_2\text{C}_2\text{O}_4 = \text{mol H}_2\text{C}_2\text{O}_4 \times 90.04 \text{ g/mol}$

**Solubility of oxalic acid (g/100 mL) =  $\frac{\text{g oxalic acid}}{\text{mL oxalic acid}} \times 100$**


# Solubility

- The maximum amount of solute that can be dissolved in a given quantity of solvent at a specific temperature.




Therefore it is the **concentration** of solute in a **saturated solution**.

- **Units:** ***g/100 mL*** (g of solute/100 mL of solution)  
***moles/L*** (moles of solute/liter of solution)

# Hydration



**Hydration** is the process in which an ion is surrounded by water molecules arranged in a specific manner.



## Heat of solution: $\Delta H_{(\text{sol})}$

- Heat of solution,  $\Delta H_{(\text{sol})}$ , is the heat generated or absorbed when a certain amount of solute dissolves in a certain amount of solvent.
- $\Delta H_{(\text{sol})}$  is positive for endothermic (heat absorbing) processes and negative for exothermic (heat generating) processes.


## Heat of solution: $\Delta H_{(\text{sol})}$ (cont'd)

$\Delta H_{(\text{sol})} = \text{lattice energy} + \text{hydration energy}$

$$\Delta H(\text{sol}) = \text{L.E} + \text{H.E}$$

- **Lattice energy (L.E):** energy **needed** to separate one mole of a solid ionic compound into gaseous ions.
- **Hydration energy (H.E):** energy **released** when an ion is surrounded by water molecules in a specific manner.

# The Solution Process for NaCl



$$\Delta H_{\text{soln}} = \text{Step 1} + \text{Step 2} = 788 - 784 = 4 \text{ kJ/mol}$$




# Conclusion

$$\Delta H(\text{sol}) = \text{L.E} + \text{H.E}$$

- **If H.E is greater than L.E**, energy is **released** and the solution process is exothermic.  
**Therefore solubility decreases with an increase in temperature.**
- **If L.E is greater than H.E**, energy is **absorbed** and the solution process is endothermic.  
**Therefore solubility increases with an increase in temperature.**

# Temperature and Solubility

## Solid solubility and temperature



solubility decreases with increasing temperature

# Experiment

- Find the solubility of oxalic acid (g/100mL) at three different temperatures.
- A saturated solution of oxalic acid at room temperature will be provided.

## a- Solubility at room temperature

- Get 20 mL of the saturated oxalic acid solution.
- Measure the temperature.
- Pipet 10 mL of the above solution into an Erlenmeyer flask and titrate against NaOH using phenolphthalein indicator

## **b- Solubility at 0° C**

### ***i - Preparation of a saturated solution at 0° C:***

- Transfer to a large test tube 20 mL of the provided saturated oxalic acid solution and 10 mL of distilled water.
- Cool the test tube in 400 mL beaker containing ice.
- Stir and wait.
- Measure the temperature (should be around 0° C).

### ***ii - Titration:***

- Pipet 10 mL of the above saturated solution carefully (leaving the solid behind) into an Erlenmeyer flask.
- Add 2 drops of phenolphthalein and titrate against NaOH.

## **c- Solubility at 40°C**

### ***i- preparation of a saturated solution at 40°C:***

- Prepare water bath at 50°C using 400 mL beaker.
- Transfer 30 mL of the sat. oxalic acid to the large test tube.
- Place the test tube in the water bath.
- Add solid oxalic acid while stirring till no more dissolves.
- Maintain the temperature around 40°C (add cold water).
- Measure the temperature.
- Decant 10 mL into a graduated cylinder, add 30 mL of water and homogenize.

### ***ii- Titration:***

- Pipet 10 mL of the above homogenized solution and titrate against NaOH using phenolphthalein indicator.
- Read the volume of NaOH and multiply by 4.