
# STRAIN GAGE TRAINING SYSTEM

### Objective



The training system offers a comprehensive introduction to the fundamentals of strain-gage technology, permitting investigation of simple mechanical situations as tension, bending and torsion.

### Fundamentals:

With metallic strain gauges, a thin metal strip or metal wire is mechanically strained inducing a change in electrical resistance. The sensitivity, i.e. ratio of change in resistance to strain, is designated by k.

$$k = \frac{\Delta R/R}{\varepsilon}$$

Where R is the resistance and  $\varepsilon$  is the strain.

In order to provide high sensitivity, one or more strain gauges are combined in a Wheatstone bridge. The bridge is supplied with a regulated DC voltage. Strain is then applied, and the occurring voltage difference (U<sub>A</sub>/U<sub>E</sub>) is displayed in a differential amplifier.

Refer to the mechanics of materials textbook, R.C Hibbler, "Mechanics of Materials", 2nd edition, Section 3.1.

#### Materials

The following material specifications shall be used in this experiment:

l: Properties of materials used in this experiment.

| contarials used in this experiment |  |          |                    |
|------------------------------------|--|----------|--------------------|
|                                    | Table 1: Properties of materials used in this experiment |          | Poisson's Ratio, µ |
|                                    |  |          |                    |
|                                    | Material   | (N/mm2)  |                    |
|                                    | Muser  | (IV/mm2) | 0.28               |
|                                    |  |          |                    |
|                                    |  | 210,000  | 0.305              |
|                                    | Steel  | 191,000  | 0.33               |
|                                    | Steel  | 191,000  |                    |
|                                    | Aluminum   | 123,000  | 0.33               |
|                                    |  | 1237     |                    |
|                                    | Copper   | 88,000   |                    |
|                                    | Brass  |          |                    |
|                                    | Blass  |          |                    |

PART 1
Set up TENSILE EXPERIMENT.




Figure 1: Sketch and photo of the setup for Group 1.

## Procedure

Ldo Gd

TECHNICAL

- Fit the tension bar in the frame as shown in figure 1. Load the bar in tension and measure the output signal to the feed voltage

- $U_a/U_e = \text{output signal to the feed voltage ratio} = \text{reading in mV/V} 10^3$ 3.



Theoretically;

Under tensile loads, stresses and strains are constant through the crosssection. Tensile stress is the ratio of the applied tensile force to the cross sectional area. The strain is related to the stress according to Hooke's law:

$$a = E^{\varepsilon}$$

- Stress = 
$$\sigma = \frac{F}{A}$$
  
-  $\sigma$  in N/mm<sup>2</sup>

F = applied load in N

A = Cross section of the bar in mm<sup>2</sup>

#### Experimentally;

Two strain gauges are attached to the opposite sides of the specimen in both longitudinal and transverse directions - as shown in figure 2. The strain gauges form the bridge shown in figure 3. Such configuration gives the following:





Figure 2: Layout of strain gauges.



Strain = 
$$\varepsilon = \frac{1}{2(1+v)} \frac{4}{k} \frac{U_a}{U_e}$$

v = Poisson's ratio- k =Sensitivity of the strain gage = 2.05




Figure 3: Bridge configuration for tensile test.

# Reference



R.C Hibbler, "Mechanics of Materials", 2nd edition, Section 3.2 and

### Keywords




Axial load, Hooke's law.

#### Report

#### 1. Discussion Questions:

- What is the Poisson's ratio? (physical properties)
- ✓ Comment on the slopes of  $\sigma$  versus  $\epsilon$  plots of different materials. Refer to graph 2 and graph 3.

#### 2. Graphs:





- ✓ Graph 2:  $\sigma$  versus  $\epsilon$  for the four bars, experimental
- ✓ Graph 3:  $\sigma$  and  $\varepsilon$  (10<sup>5</sup>) versus material, column chart, experimental, under 50N load.
- ✓ Comment on Hooke's law.



PART 2
Set up

BENDING TEST.





### Procedure

Les Gel TECHNICAL

Fit the tension bar in the frame, connect up, switch, and adjust the

2. Load the bar in bending and measure the output signal to the feed voltage rano.  $U_{\text{B}}/U_{\text{e}} = \text{output signal to the feed voltage ratio} = \text{reading in mV/V } 10^{-3}$ 

# Theoretically.

Under benefity, the beam publishes maximum strongs a the surface.

o m himm


M = Moment due to applies tout a % am-

I - Moment of merita of the Cross section of he has a new  $\xi = Largest$  distance from the sentral axis to the surface


#### Experimentally:

Two strain gauges are provided in each of the amount and compression sides. Unlike the tension test, the four strain gauges are arranged a fielongitudinal direction - figure 5. The resulting widge a down in Squee 8.





& = Sensitivity of the serum page = 2.85



# Reference



R.C Hibbler, "Mechanics of Materials", 2nd edition, from Section 6.3 &

### Keywords



Bending, Neutral axis, Flexure formula.

#### Report

- 1. Discussion Questions:
- ✓ What are c and I in the stress expression? Explain.
- ✓ Analyze the effect of load position on stress.




✓ σ versus x for steel beam, under 20N load, experimental and theoretical.

PART 3

Set up

TORSION TEST.






Figure 7: Sketch and photo of the setup for Group 3.

Procedure

1. Load the bar in torsion and measure the output signal to the feed voltage ratio = reading in mV/V 10<sup>3</sup>

Voltage ratio.  $U_8/U_e = \text{output signal to the feed voltage ratio} = \text{reading in mV/V 10}$ 

Theoretically;

The torsional moment is the product of the applied force and the torsional

- F = Applied load (N)
- Arm = torsional arm (mm)



The torsion bar is provided with four strain gauges located along the direction of the principal normal stresses, i.e. at 45° with the neutral axis.




Figure 8: Layout of strain gauges, Bar under torsion.

The torque in the bar is calculated from the shear stress measured as

- As noted previously, the strain is  $\varepsilon = \frac{1}{k} \frac{U_a}{U_c}$ follows:
- Assuming pure shear stress, the shear strain is twice the normal strain:
- According to Hooke's law, shear stress and shear strain are related as
- follows:  $\tau = G.\gamma$
- Torsional Moment (Torque) =  $M_t = \tau W_p$ 
  - Therefore  $M_1 = G \frac{2}{k} \times \frac{U_a}{U_s} \times W_p$ k =Sensitivity of the strain gage = 2.05
- $W_p =$ Section modulus of torsion for the circular cross section G =Shear modulus of steel =  $80,000 \text{ N/mm}^2$

# Reference



R.C Hibbler, "Mechanics of Materials", 2nd edition, from Section 5.1 &

# Keywords



Torque, Shear stress, Shear strain, Shear modulus, Section Modulus The torsion formula.

#### Report





- ✓ Experimentally, how the load (F) affects the torsional moment (M₁).

  Assess with the theory.
- 2. Graph:
- ✓ M<sub>t</sub> versus F for steel beam, experimental and theoretical.